atect

Date of issue 2026/9/25

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atect

Date of issue 2026/9/25

atect Stackable Medium for coliform count “aS-Medium CC” instruction manual

Product Information

aS-Medium CC consists of a plastic substrate and a cover film with a resin frame. This product is an easy-to-use microbiological test medium that measures the coliform counts from a sample solution as an alternative to conventional agar-based methods. The plates can be stacked stably due to the indentations at the four corners. By adding an adhesive resin frame to the cover film, the product improves adhesion between the cover film and the substrate, prevents the cover film from unintentionally turning over, and prevents leakage of the test liquid. In addition, because the product is in the form of a thin sheet, it contributes to space saving during storage and reduces environmental impact during disposal.

aS-Medium CC is intended for use by trained laboratory specialists or by general users under the supervision of such specialists in a microbiology laboratory setting.

PTM Certification

aS-Medium CC has been evaluated against a reference method (BAM Ch. 4) in the AOAC Research Institute PTM program , and is currently being prepared for AOAC PTM certification application for counting coliform bacteria in meat products (ready-to-eat ham, ready-to-eat bacon, fresh raw pork, fresh raw beef trim, and fresh raw chicken). For evaluation in this test, 2 sheets are used for each test area, and the number of bacteria is calculated based on the average value and compared with the reference method.

Materials Required, not Provided

1. Incubator (35±1℃).
2. Stomacher.
3. Stomacher bag with filter (250 mL or 500 mL volume).
4. Butterfield’s phosphate-buffered dilution water.
5. Pipettor, pipette, and tips.

Sample Preparation

1. Weigh 50 g of the sample and add it to the filtered Stomacher bag.
2. Add 450mL of the appropriate buffer solution (such as Butterfield's phosphate-buffered diluent) to the sample to filtered Stomacher bag.
3. Homogenize the sample using a Stomacher.
 Note: Alternatively, for steps 1 to 3 above, weigh two 25 g portions into separate Stomacher bags. Add 225 mL of Butterfield's phosphate-buffered diluent to each portion and homogenize for 2 minutes. After homogenization, combine the two portions into a single container and briefly mix by hand.
4. Dilute the sample solution as needed.

Inoculation

1. Pinch the cover film between your thumb and index finger and place your remaining fingers from the center of the film over the logo in the upper left corner.
2. Using your thumb and index finger, lift and turn over the bottom edge of the resin frame of the cover film.
3. While holding down the product, turn over the cover film.
4. Using a pipette, dispense 1 mL of the sample solution onto the medium. Slowly replace the cover film to cover the surface of the medium, taking care not to trap air bubbles. Let stand for 3 minutes. (If the product is moved or tilted before gelation occurs, the sample solution may leak out of the culture layer.)
5. Stack the plates with the cover film on top. Up to 25 plates may be stacked.
6. Incubate the plate at 35±1℃ for 24±2h.
7. If necessary, the cover film can be peeled off to perform inoculation.

Interpretation

The medium contains the enzyme substrate X-Gal (5-Bromo-4-Chloro-3-Indolyl-β-D-Galactoside). The coliform bacteria colonies turn blue. The number of coliform bacteria can be measured by counting the blue colonies.

Conduct the test using two sheets per test area in accordance with the test method used for PTM certification. Calculate the average number of bacteria per plate, and determine the number of bacteria in the sample based on this value.

Precautions for counting colonies

The area of the medium is about 20 cm^2. The range of the number of bacteria that can be counted is 30 to 300 cfu. If the number of colonies is large, it can be estimated by counting 1 square cm (1 cm x 1 cm) in the grid and multiplying by 20.

Storage and Shelf life

1. Store the product in a refrigerator at 4–10 °C.
2. The shelf life is 12 months when stored at 4–10 °C and kept unopened.
3. The expiration date is printed on the lower left of the product, as well as on the label on the aluminum bag.
4. After opening, return the product to the aluminum bag, seal it, and store it in a refrigerator at 4–10 °C. Use it within 1 month.

Precautions for use

1. Check the instruction manual before using this product.
2. This product is mainly used for microbiological testing of food and environmental samples. Do not use it for purposes other than its intended use.
3. The product should only be handled by experienced individuals.
4. This product has not been verified for all foods, inspection procedures, and strains.
5. This product detects coliform bacteria using β-galactosidase as an indicator. Some non-coliform bacteria (such as Serratia species) may also be detected as coliform bacteria.
6. If a food sample contains components that react with the chromogenic enzyme substrates, the entire culture sheet may become discolored. In such cases, further dilute the sample before inoculation.
7. Please note that the user is responsible for the interpretation and operation of the inspection results, and the manufacturer assumes no responsibility for these activities.
8. Products that have passed the expiration date are not covered by the warranty.
9. If the product is damaged, contains foreign matter, or is discolored, refrain from using it and contact the manufacturer or distributor.

Safety Precautions

1. Assay users should observe standard microbiological practices and safety precautions when performing this assay.
2. If the product or reagents come into contact with the eyes or mouth, rinse thoroughly with water and consult a physician.
3. Used media must be disposed of after autoclaving. It is the responsibility of the end user to handle waste and effluents according to their nature and degree of hazard, and to treat and dispose of them in accordance with applicable local, state, and federal regulations.

Manufacturer information

Atect Corp.
Address: 3275-1, Kamihaneda-cho, Higashiomi city, Shiga,
527-0082, Japan
TEL: 0748-20-3535 FAX: 0748-25-7171
Mail: info-bio@atect.co.jp